













# PURESPIN PRP™

**\*\* PLEASE DISCARD RED VENTED CAP FROM CONCENTRATING DEVICE BEFORE USE \*\***

Note: Always swab self-sealing port with sterile alcohol prior to accessing with a sterile syringe

<p><b>Step 1:</b></p>  <p>Draw 5mL of Citrate Anticoagulant into 30mL Syringe</p>	<p><b>Step 2:</b></p>  <p>Draw 25mL whole blood from the patient, filling the syringe to 30mL</p>	<p><b>Step 3:</b></p>  <p>Load anticoagulated whole blood into the <b>Concentrating Device</b></p>	<p><b>Step 4:</b></p>  <p>Counterbalance and process the <b>Concentrating Device</b> at</p> <p><b>1 minute 4200 RPM</b></p>
<p><b>Step 5:</b></p>  <p>Using the 30mL syringe, aspirate the platelet plasma suspension (PPS) until RBC fills the aspirating pipe.</p> <p>(Its normal to aspirate small amounts of RBC into the syringe during this process)</p>	<p><b>Step 6:</b></p>  <p>Transfer the platelet plasma suspension (PPS) into the <b>Concentrating Accessory</b></p>	<p><b>Step 7:</b></p>  <p>Counterbalance and process the <b>Concentrating Device</b> at</p> <p><b>5 minutes 4200 RPM</b></p>	<p><b>Step 8:</b></p>  <p>Platelet Concentrate Buffycoat</p> <p>Platelet concentrate buffycoat separates out at the bottom of the <b>Concentrating Accessory</b></p>
<p><b>Step 9:</b></p>  <p>Aspirate Platelet Poor Plasma</p> <p>Leave approximately 4mL</p> <p>Aspirate platelet poor plasma from the <b>Concentrating Accessory</b> Leave 4mL of plasma.</p>	<p><b>Step 10:</b></p>  <p>Attach the 12mL syringe and swirl to resuspend the platelet buffycoat into the plasma.</p>	<p><b>Step 11:</b></p>  <p>Tilt to immerse Aspirating Pipe into Purespin PRP®</p> <p>Tilt to immerse the Aspirating Pipe into the Purespin PRP®</p>	<p><b>Step 12:</b></p>  <p>4mL Purespin PRP®</p> <p>Extract the Purespin PRP® into the 12mL syringe.</p>